
Multiplanter
Cartridge

Multiplanter

New world of surgery unfolds


INRUT

GNU Multiplanter Introduction

A secret method of
lifting only once and planting ten times

Inventor

Structural design development Medical Team, Ohdae Research Institute

How did you come to develop a Multi-hair implanter?

Ohdae Research Institute has been researching robotic and automatic hair transplant device for many years. In the course of our research, we invented a device that rotates automatically without electric power and significantly improves the speed of the surgical procedure. We have developed it to be as simple as possible and highly efficient so that users can quickly adapt to using it.


Developed a clinical surgical method <Daegu> 'K' university hospital Hair Transplantation Center Professor Kim Moon-kyu

What are the advantages of the surgical method using the multi-hair implanter?

The multi-hair implanter is faster than anything else. The multi-hair implanter can be loaded with ten hair follicles at a time allowing continuous transplant. Therefore, the implant time is reduced by more than an hour, and the procedure is that much less painful for the patients. Since the implant time is reduced, we are expecting improved engraftment rate from shortened time of hair follicles exposed in vitro.

*President of the Korean Society of
Hair Restoration Surgery

NEW WORLD OF SURGERY UNFOLDS

How do you feel about using the multi-hair transplanter?

*Vice President of the
Korean Society of Hair Restoration Surgery

Director Ahn Ji-Sup Q&A

(<Seoul> 'A' clinic)

In addition to reducing the procedure time, INRUT also raised the satisfaction of assistants (loaders) participating in the procedure. It has become more relaxing for the assistants during the procedure because they can load ten hair follicles once they hold an implanter. In the past, the procedure took a long time and often required the assistants to take turns during the procedure. However, thanks to reduced procedure time with, we could put an end to assistant rotation and enhance our concentration level during the procedure.

Director Jung Jin-wook Q&A

(<Gwangju> 'M' clinic)

Compared to existing hair implanters, using the multi-hair transplanter significantly reduced the procedure time. The treatment can be performed within two hours for 3,000 hair follicles, and it only takes about two and half hours for 4,000 hair follicles. With shorter procedure time, patients go through less pain, which can increase their satisfaction with the procedure. Nurses who are responsible for loading hair follicles prefer this multi-hair implanter. I think this is a great product that is beneficial to practitioners, assistant nurses, and patients.

* This product may be used only by authorized medical personnel. / This contents was made only for medical personnel. / This contents include a user's personal opinion of the product.

Q&A

What about the hair follicle loading time?

With the existing hair transplanters, we had to load one hair follicle on a device at a time. But the multi-hair implanter comes with a cartridge that can be rotated and loaded with ten hair follicles quickly, which significantly reduced the loading time.

Q&A

How different from conventional hair transplant procedure in terms of speed?

The transplanting time, which previously took about two and half hours, has been reduced to about one and half hour. We used to have two or three people do the loading and one person pass the device to reduce the procedure time. But with the multi-hair implanter, I think we can further reduce the time with ust two people loading.

Q&A

What about its conveniences such as use of anesthetics and assistant rotation time?

The use of anesthetics has also been reduced than before, and the number of loading assistants and delivery assistants has also decreased. Currently, we only need two loading assistants. When using a single-hair transplanter, I felt great fatigue from the repetitive task of planting a hair and replacing the device one after another. But it is much more conven ient now since I can implant ten follicles in a row with one replacement.

Director Park Geun Q&A

(<Busan> 'S' clinic)


I was pleased with the fact that INRUT can reduce anesthesia time and the amount of anesthetic used on patients. I have to administer additional anesthetic in the middle of the procedure because the patients complain of pain when the anesthesia administered at the beginning of the procedure wears off after a period of time. With the use of the multi-hair transplanter, the procedure time was reduced, and so did the anesthesia time, the number of administering anesthetics and amount of anesthetics used. Consequently, INRUT helped relieve thereby my patients' physical burden.

Professor Lee Yong-hyun Q&A

(<Daegu> 'K' university hospital)

Other people who use a multi-hair transplanter often use three units, but I have no problem conducting the procedure with two units. This is because the loading speed of my assistant who is in charge of loading hair follicles has gotten faster compared to the loading speed with a general hair implanter. When conducting the transplant procedure with two units of INRUT, the number of assistants is reduced by one, which can increase the efficiency of the procedure.


5-step hair transplant method using INRUT multi-hair transplanter


STEP. 1

Grabbing the device
Hold the handle with your thumb and middle finger, and gently place your index finger on the rear cap.


STEP. 2

Stick needle
Insert needle into the implant spot.

STEP. 3


Remove needle
Lift the sliding handle to remove needle out of the skin. At this time, it is advisable to lift the handle slightly before it clicks.

STEP. 4

Remove nozzle from scalp
Lift up the hair transplanter to remove the nozzle from the skin. Check the condition of the implanted hair follicle.

(*Arranging the hair follicles is recommended after all 10 hair follicles have been implanted.)

STEP. 5

Rotate cartridge
Raising the sliding handle all the way up allows the cartridge to rotate. Loosen your index finger completely to lower the sliding handle all the way down.
(*Take note that the cartridge will not rotate if the sliding handle is not raised all the way up until it clicks.)
(*If the sliding handle is not lowered completely, the metal pin rod will not get caught in the latch, causing an error during implant.)

hair transplant method using
INRUT multi-hair transplanter
QR CODE

INRUT ERROR CODE

ERROR CODE. 1

When inserting the needle into the skin, be careful not to press the rear cover with your index finger.

ERROR CODE. 2

If the nozzle slips off the skin due to a latch release jolt when removing the needle out of the skin, the nozzle must be lifted from the scalp before the sliding handle clicks.

ERROR CODE. 3

When the sliding handle is lowered again after implanting a hair follicle, the next needle sticks out of the nozzle. Be careful not to accidentally pierce the scalp with the needle.

ERROR CODE. 4

If you do not lift the sliding handle all the way up, the cartridge will not rotate. The handle must be lifted all the way up until it clicks. Also, it must be lowered all the way down when lowering.

※ How to check the normal position of stopper ring: Check whether the metal pin rod temporarily stays at the tip of the nozzle when the sliding handle is slowly raised.

ERROR CODE. 5

If the stopper ring is placed too high, the metal pin rod will not work properly. It will not temporarily stay at the tip of the nozzle when the handle is raised. Therefore, make sure to place the stopper ring at the correct position according to the adjustment method.

ERROR CODE. 6

If the stopper ring is placed too low, the metal pin rod will stick out too much from the nozzle. This will affect the depth of implant, so make sure to place the stopper ring at the correct position according to the adjustment method.


ERROR CODE. 7

The slit channel on the needle comes in two types: normal and wide. In general, the normal type is appropriate for hair follicle that is trimmed and arranged at the upper part, while the wide type is suitable for untrimmed hair follicle. Select a channel that matches the hair follicle shape so that the hair follicle does not slip off the channel during implant.


ERROR CODE
QR CODE

How to replace the cartridge in INRUT multi-hair transplanter

**STEP. 1****Remove nozzle**

Unwind the stopper ring to remove it. Then, pull back the sliding handle and open the case joint using a snap pliers to remove the nozzle.

**STEP. 2****Combine cartridge cover**

Pull back the handle and close the cartridge cover.

**STEP. 3****Remove rear cover, pusher, and cartridge**

Open the latches on the rear cover one by one to remove the rear cover. Then, remove the springs and the pusher, and keep them aside.

**STEP. 6****Combine pusher, spring, and rear cover**

Insert the pusher and two springs into the case aligning them in the right position, and then close the rear cover. Then, fasten the latches and check whether the rear cover is properly combined.

**STEP. 7****Remove cartridge cover**


Twist open the cartridge to remove it. (*Be careful not to injure your hand on the sharp edge of the case.)

**STEP. 8****Combine nozzle**

With the sliding handle raised, combine the new nozzle to the case.

* Caution: Use the case as a base tray to prevent loss of parts.


Search INRUT on YouTube and find out more.

**STEP. 4****Clean and sterilize parts**


Remove the used cartridge, and then clean and sterilize the parts that are to be reassembled. (*Use sterilizing method for non-critical items)

**STEP. 5****Prepare new cartridge and nozzle**


Prepare a new needle cartridge and a nozzle. Insert the new needle cartridge into the case.

**STEP. 9****Combine stopper ring**

Wind the stopper ring all the way in to combine it with the case, and then loosen it two notch in the opposite direction to place it in the correct position. (*Malfunction may occur if the stopper ring is not placed in the correct position.)

**STEP. 10****Check metal pin rod in nozzle**

(*When removing the nozzle cap, pull back the sliding handle to prevent the damage of needles.) With the nozzle cap taken off, slightly raise the sliding handle to check whether the metal pin rod stays at the tip of the nozzle.


**Cartridge replacement video****Scan QR code**

You can check the video of replacing the INRUT multi-hair transplanter cartridge on your mobile device.

INRUT needle type

* Can be custom made by selecting type of needle according to intended use

Ways to sterilize INRUT multi-Implanter

- The manufacturer delivers all products in a sterilized condition after EO gas processing at the time of shipment.
- The manufacturer of this device specifies the EO gas sterilization method as the sterilization method for this product.
- Do not expose this product to high temperatures. [This product may not be sterilized in an autoclave.]
- Reusable instruments and items should be thoroughly cleaned before disinfection or sterilization.
- Guidelines for other disinfection and sterilization methods are as follows. [Guidelines by Korea Centers for Disease Control and Prevention]

Needle cartridge - Critical items

- E.O.(Ethylene oxide) gas
- Hydrogen peroxide gas plasma
- 2% or higher concentration glutaraldehyde [Cidex, Glutaral, Widex, Totacid-28] 20~25°C, 10 hours
- 7.5% Hydrogen peroxide, 6 hours
- When using peracetic acid, be extra cautions about metal corrosion.

Nozzle - Semicritical items

- Glutaraldehyde [Cidex, Glutaral, Widex, Totacid-28] 2%, 20°C, 20 minutes / 2.5%, 35°C, 5 minutes
- 0.55% Ortho-phthalaldehyde (OPA) [Cidex-OPA] 20°C, 12 minutes / 25°C, 5 minutes
- 7.5% Hydrogen peroxide 30 minutes
- 650~675 ppm or higher concentration Hypochlorite, 10 minutes [Product made by electrolysis at the place of use]

Other parts of case - Non-critical items

- Sodium hypochlorite with effective chlorine content of 100 ppm or more for 1 minute or longer
- Phenolic disinfectants for 1 minute or longer
- Iodophor disinfectants for 1 minute or longer
- Quaternary Ammonium Chloride for 1 minute or longer
- 70~90% alcohol (Ethanol/Isopropanol) for 1 minute or longer

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